

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Formaldehyde

Parameter*	Value*
REF	91328
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, additional reagents, instruction leaflet in a box
Parameter	Formaldehyde
Gradation	0 • 10 • 20 • 40 • 60 • 100 • 200 mg/L HCHO
Method	Triazol
Color change	beige → blue-violet
Instrumental measuring range	10–200 mg/L HCHO

*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG
Valenciener Str. 11
52355 Düren · Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
US Tel.: +1 888 321 62 24 sales-us@mn-net.com