

## Data sheet

## Filter paper

| Parameter*                    | Value*                                                                      |
|-------------------------------|-----------------------------------------------------------------------------|
| Grade                         | MN 618                                                                      |
| Basis weight                  | 140 g/m <sup>2</sup>                                                        |
| Thickness                     | 0.33 mm / 13.0 mil                                                          |
| Filtration speed              | 23 s                                                                        |
| Surface                       | Smooth                                                                      |
| Color                         | White                                                                       |
| Material                      | Cellulose                                                                   |
| Ash content                   | ca. 0.1 %                                                                   |
| Characteristics / Application | Thick filter paper, medium fast filtration, for large quantities of liquids |

\*These data are typical values and must not be regarded as specification.