

Receiver Plates

1. Contents

Receiver Plates 10 µm	
REF	740989.4
Receiver Plates 10 µm	4
Receiver Plates 35 µm	
REF	740512.4
Receiver Plates 35 µm	4
Receiver Plates 35 µm hydrophilized	
REF	740513.4
Receiver Plates 35 µm hydrophilized	4
Receiver Plates 50 µm	
REF	740688.4
Receiver Plates 50 µm	4
Receiver Plates 50 µm hydrophilized	
REF	740689.4
Receiver Plates 50 µm hydrophilized	4

2. Product description

Receiver Plates are 96-well plates available with inserted filter frits of 10 µm, 35 µm, or 50 µm pore size. They can be used for general filtration purposes as well as for retaining chromatographic resins (e.g., Protino® Ni-NTA Agarose, Protino® Glutathione Agarose 4B, Sephadex® G25, G50, or Sephacryl® S200). The Receiver Plates can be used under centrifugation using standard benchtop centrifuges as well as under vacuum. The plates are also offered with hydrophilized filter frits which enable the plates to be used additionally with gravity flow.

3. Specifications

Plate material:	Polypropylene
Filter frit material:	Polyethylene untreated (for use under centrifugation and vacuum) Polyethylene hydrophilized (for use under gravity, centrifugation, and vacuum)
Filter frit pore size:	approx. 10 µm, 35 µm, or 50 µm
Filter frit diameter:	8.3 mm
Effective filtration area:	54.1 mm ²
Volume capacity per well:	1.5 mL

For research use only.

4. Safety instructions

Receiver Plates do not contain hazardous contents.

5. Product use restriction / warranty

All MACHEREY-NAGEL products are designed for their intended use only. They are not intended to be used for any other purpose. The description of the intended use of the products can be found in the original MACHEREY-NAGEL product leaflets. Before using our products, please observe the instructions for use and the safety



instructions from the respective Material Safety Data Sheet of the product.

This MACHEREY-NAGEL product is carrying documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. The provided warranty is limited to the data specifications and descriptions as given in the original MACHEREY-NAGEL literature. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agents or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized. They should not be relied upon by the customer and are not a part of a contract of sale or of this warranty.

Liability for all possible damages that occur in any connection with our products is limited to the utmost minimum as stated in the general business terms and conditions of MACHEREY-NAGEL in their latest edition which can be taken from the company's website. MACHEREY-NAGEL does not assume any further warranty.

Products and their application are subject to change. Therefore, please contact our Technical Service Team for the latest information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Descriptions in MACHEREY-NAGEL literature are provided for informational purposes only.

Last updated: 08/2022, Rev. 04

Please contact:

MACHEREY-NAGEL GmbH & Co. KG

Tel.: +49 24 21 969-333

support@mn-net.com

Trademarks:

Protino is a registered trademark of MACHEREY-NAGEL GmbH & Co. KG

Sephadex and Sephacryl are registered trademarks of GE Healthcare companies