


NANOCOLOR® Organische Komplexbildner 10

Organic Complexing Agents 10 / Complexants organiques 10



Test 0-52
REF 985052







Methode / Method / Méthode / Método

560 / 540 nm

0521 0.5 – 15.0 mg/L I_{BiK}
0522 0.5 – 20.0 mg/L EDTA


Min	20 – 80 %		Max
0.5	3.4	12.1	15.0 mg/L
0.5	4.4	16.1	20.0 mg/L

Nullwert erforderlich. Siehe Rückseite. / Blank value required. See back. / Blanc. Voir au verso. / Valor en Blanco. Ver reverso.



NANOCOLOR® Organische Komplexbildner 10

Organic Complexing Agents 10/Complexants organiques 10



Test 0-52
REF 985052

Nullwert/Blank value/Blanc/Valor en Blanco

