


visocolor® ECO Eisen 2


Iron / Fer / Hierro

Messbereich / Range / Domaine de mesure / Rango

0.04 – 1.0 mg/L Fe




Probe
Sample
Echantillon
Muestra




4 ∅ Fe-1




Schütteln
Shake
Agiter
Agitar




1 ∅ Fe-2



Schütteln
Shake
Agiter
Agitar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser

Suitable for seawater

Convient à l'eau de mer

Adecuado para el agua de mar