



# visocolor® ECO Kupfer

Copper / Cuivre / Cobre

Messbereich / Range / Domaine de mesure / Rango  
0.1 – 1.5 mg/L Cu<sup>2+</sup>




1  
Probe  
Sample  
Echantillon  
Muestra




2  
5 ∅ Cu-1




3  
Schütteln  
Shake  
Agiter  
Agitar




4  
5 ∅ Cu-2



5  
Schütteln  
Shake  
Agiter  
Agitar



6  
Warten  
Wait  
Attendre  
Espere



7  
Messung  
Measurement  
Mesure  
Medición

Geeignet für Meerwasser  
Suitable for seawater  
Convient à l'eau de mer  
Adecuado para el agua de mar

