



# visocolor® ECO Cyanid

Cyanide / Cyanure / Cianuro

Messbereich / Range / Domaine de mesure / Rango  
0.01 – 0.20 mg/L CN<sup>-</sup>




1  
Probe  
Sample  
Echantillon  
Muestra




2  
5 ∅ CN-1




3  
Schütteln  
Shake  
Agiter  
Agitar



4  
1 ∅ CN-2




5  
Schütteln  
Shake  
Agiter  
Agitar




6  
5 ∅ CN-3



7  
Schütteln  
Shake  
Agiter  
Agitar



8  
Warten  
Wait  
Attendre  
Espere



9  
Messung  
Measurement  
Mesure  
Medición

Geeignet für Meerwasser nach Verdünnung 1 + 3  
Suitable for sea water after dilution 1 + 3  
Convient à l'eau de mer après dilution 1 + 3  
Adecuado para el agua de mar tras su dilución 1 + 3



MACHEREY-NAGEL  
www.mn-net.com