

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Phosphate

Parameter*	Value*
REF	91320
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, additional reagents, instruction leaflet in a box
Parameter	Phosphate
Gradation	0 • 3 • 10 • 25 • 50 • 100 mg/L PO ₄ ³⁻
Method	Phosphomolybdenum blue
Color change	white → blue-green
Instrumental measuring range	3–80 mg/L PO ₄ ³⁻

*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG
 Valencienner Str. 11
 52355 Düren · Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com