

Data sheet

Filter paper

Parameter*	Value*
Grade	MN GF-3
Basis weight	52 g/m ²
Thickness	0.27 mm / 10.6 mil
Filtration speed	25 s
Maximum temperature	500 °C / 932 °F
Surface	Smooth
Color	White
Material	Borosilicate glass micro fibers
Average retention capacity	0.6 µm
Characteristics / Application	Glass fiber filter without binder, e.g. for laboratory filtration and air sampling

*These data are typical values and must not be regarded as specification.