



visocolor® ECO Cyanid

Cyanide / Cyanure / Cianuro

Messbereich / Range / Domaine de mesure / Rango
0.01 – 0.20 mg/L CN⁻




Probe
Sample
Echantillon
Muestra




5 ∆ CN-1




Schütteln
Shake
Agiter
Agitar




1 ∆ CN-2




Schütteln
Shake
Agiter
Agitar




5 ∆ CN-3



Schütteln
Shake
Agiter
Agitar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser nach Verdünnung 1 + 3
Suitable for sea water after dilution 1 + 3
Convient à l'eau de mer après dilution 1 + 3
Adecuado para el agua de mar tras su dilución 1 + 3