


Data sheet

Filter paper

| Parameter* | Value* |
|-------------------------------|---|
| Grade | MN 616 LSA-50 |
| Maximum temperature | 100 °C / 212 °F |
| Surface | Smooth |
| Color | White |
| Material | Cellulose, cation exchanger |
| Capacity | 0,05 ± 0,01 meq/cm ² , 0,9 ± 0,2 meq/Filter with ø 48 mm |
| Characteristics / Application | Filter paper with strongly acidic cation exchange resin, Matrix polystyrene crosslinked with DVB, functional groups SO ₃ H, supplied in H ⁺ form/ |

*These data are typical values and must not be regarded as specification.


 18.02.2026 · Rev. 1 · www.mn-net.com

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com