


Data sheet

Filter paper

| Parameter* | Value* |
|-------------------------------|---|
| Grade | MN 210 PE |
| Basis weight | 140 g/m ² |
| Thickness | 0.25 mm / 9.8 mil |
| Surface | Smooth |
| Color | White |
| Material | Cellulose, polyethylene |
| Characteristics / Application | Filter paper, one side polyethylene (PE)-coated, e.g. for covering laboratory workbenches |

*These data are typical values and must not be regarded as specification.



18.02.2026 · Rev. 1 · www.mn-net.com

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com
CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com
US +1 888 321 62 24 sales-us@mn-net.com