

MACHEREY-NAGEL

# Thin layer chromatography summer promotion

30 %  
off\*

Chromatography



State-of-the-art TLC products

- Efficient
- Reproducible
- Robust

Timer with  
purchase



**MACHEREY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)



# TLC summer promotion



## What is TLC?

It's a quick, easy and cost-efficient method for separating compound mixtures utilizing planar plates or sheets made of glass, aluminum or polyester coated with silica or other separating media.

## What is it used for?

A common method in pharmacy, food industries, natural compound analyses or organic syntheses for control of identity, yield or composition of substances and products.

## Why should customers buy MACHEREY-NAGEL TLC products?

- Comprehensive range of plate sizes, surface chemistries and backings
- Pre-coated plates ready for immediate use
- Homogeneous, smooth and well adhering layers
- Consistent high quality from batch-to-batch and from plate-to-plate

Timer with purchase



## Ordering information

Designation	REF	Plate size (cm)	Thickness of layer (mm)	Pack of	List price per pack (\$)	Promotion price per pack (\$)
<b>ALUGRAM® Xtra aluminum sheets</b>						
■ ALUGRAM® Xtra SIL G/UV <sub>254</sub>	818329	2.5 x 7.5	0.2	200	<del>206.00</del>	144.20
■ ALUGRAM® Xtra SIL G/UV <sub>254</sub>	818331	4 x 8	0.2	50	<del>55.40</del>	38.78
■ ALUGRAM® Xtra SIL G/UV <sub>254</sub>	818332	5 x 20	0.2	50	<del>112.00</del>	78.40
■ ALUGRAM® Xtra SIL G/UV <sub>254</sub>	818333	20 x 20	0.2	25	<del>203.00</del>	142.10
<b>Glass plates</b>						
■ TLC plates SIL HD UV <sub>254</sub>	809220	10 x 10	0.25	25	<del>116.00</del>	81.20
■ TLC plates SIL HD UV <sub>254</sub>	809222	10 x 20	0.25	50	<del>305.00</del>	213.50
■ TLC plates SIL HD UV <sub>254</sub>	809223	20 x 20	0.25	25	<del>270.00</del>	189.00
■ TLC plates SIL HD UV <sub>254</sub>	809227	5 x 10	0.25	50	<del>208.00</del>	145.60
<b>POLYGRAM® polyester sheets</b>						
■ POLYGRAM® SIL G/UV <sub>254</sub>	805901	2.5 x 7.5	0.2	200	<del>206.00</del>	144.20
■ POLYGRAM® SIL G/UV <sub>254</sub>	805021	4 x 8	0.2	50	<del>55.40</del>	38.78
■ POLYGRAM® SIL G/UV <sub>254</sub>	805022	5 x 20	0.2	50	<del>112.00</del>	78.40
■ POLYGRAM® SIL G/UV <sub>254</sub>	805023	20 x 20	0.2	25	<del>203.00</del>	142.10

Promo code: **TLC26**

Quote **TLC26** when ordering by e-mail, phone, fax or through our webshop.

\* Special promotion prices only valid for the mentioned products.

Offer valid in the USA only from July 1<sup>st</sup> until September 31<sup>st</sup> 2026. Quote TLC26 when ordering. No other discounts apply.



[www.mn-net.com](http://www.mn-net.com)

# MACHEREY-NAGEL



MACHEREY-NAGEL Inc. · 924 Marcon Blvd., Suite 102 · Allentown, PA 18109 · USA

US +1 888 321 62 24 sales-us@mn-net.com  
DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com  
FR +33 388 68 22 68 sales-fr@mn-net.com