

Data sheet

31.10.2022

Filter paper

Parameter*	Value*
Grade	MN 616 LSB
Surface	smooth
Maximum temperature	60 °C
Capacity	0,03 ± 0,005 meq/cm ² , 0,55 ± 0,1 meq/Filter with ø 48 mm
Characteristics / Application	Filter paper with strongly basic anion exchange resin, Matrix polystyrene crosslinked with DVB, functional groups quarternary ammonium compounds, supplied in OH- form/

*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG
 Valencienner Str. 11
 52355 Düren · Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com