



visocolor® ECO DEHA

DEHA / DEHA / DEHA


Messbereich / Range / Domaine de mesure / Rango
0.01 – 0.30 mg/L DEHA




filtrieren wenn trüb
filter if cloudy
filtrer si trouble
filtra si turbido




Probe
Sample
Echantillon
Muestra



4 ∅ DEHA-1




Schütteln
Shake
Agiter
Agitar




4 ∅ DEHA-2



Schütteln
Shake
Agiter
Agitar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar

