

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Peracetic acid 500

Parameter*	Value*
REF	91341
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Peracetic acid
Gradation	0 • 50 • 100 • 200 • 300 • 400 • 500 mg/L Peracetic acid
Method	Redox reaction
Color change	yellow → green
Instrumental measuring range	50–500 mg/L Peracetic acid

*These data are typical values and must not be regarded as specification.

