





visicolor® ECO Kalium

Potassium / Potassium / Potasio


Messbereich / Range / Domaine de mesure / Rango
2 – 15 mg/L K⁺




1
Spülen
Rinse
Rincer
Lavar




2
Probe bis Markierung
Sample up to mark
Échantillon jusqu'au repère
Muestra para marcar




3
15 Δ K-1




4
Umschwenken
Shake gently
Secouer légèrement
Mezclar volteándolo




5
1 Δ K-2



6
Umschwenken
Shake gently
Secouer légèrement
Mezclar volteándolo



7
Messröhrchen füllen
Shake gently
Filling measuring tube
explétif la cuve de mesure
echar en tubito de medida



8
bis Kreuz unsichtbar
till cross invisible
jusque croix invisible
hasta cruz invisible

Geeignet für Meerwasser nach Verdünnung 1 + 1
Suitable for sea water after dilution 1 + 1
Convient à l'eau de mer après dilution 1 + 1
Adecuado para el agua de mar tras su dilución 1 + 1

