

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Glutaraldehyde

Parameter*	Value*
REF	91343
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Glutaraldehyde
Gradation	0•0.5•1.0•1.5•2.0•2.5 % Glutaraldehyde
Method	Mixed indicator
Color change	bright orange → magenta

*These data are typical values and must not be regarded as specification.

