

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Total iron 1000

Parameter*	Value*
REF	91330
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Iron
Gradation	0 • 5 • 20 • 50 • 100 • 250 • 500 • 1000 mg/L Fe ²⁺ 3 ⁺
Method	2,2'-Bipyridine
Color change	white → dark red

*These data are typical values and must not be regarded as specification.

