



visocolor® ECO Aluminium

Aluminum / Aluminium / Alumino

Messbereich / Range / Domaine de mesure / Rango
0.10 – 0.50 mg/L Al³⁺




1
Probe
Sample
Echantillon
Muestra




2
1 — Al-1




3
Schütteln
Shake
Agiter
Agitar



4
4 ∅ Al-2




5
Schütteln
Shake
Agiter
Agitar




6
10 ∅ Al-3



7
Schütteln
Shake
Agiter
Agitar



8
Warten
Wait
Attendre
Espere



9
Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser nach Verdünnung 1 + 9
Suitable for sea water after dilution 1 + 9
Convient à l'eau de mer après dilution 1 + 9
Adecuado para el agua de mar tras su dilución 1 + 9



MACHEREY-NAGEL
www.mn-net.com