Chromatography

MACHEREY-NAGEL

CHROMAFIL® Xtra Winter promotion

Buy 3 packs and get free HARIBO® gummy bears





For reliable and successful filtrations

- PTFE and GF/PTFE syringe filters
- PES syringe filters
- PA and GF/PA syringe filters
- CA syringe filters

As low as $ 0.50/filter*

CHROMAFIL® Xtra Winter promotion

Buy 3 packs and get free HARIBO® gummy bears


As low as $ 0.50/filter*




How to select the optimal CHROMAFIL® syringe filter

1. Filter size

Sample volume


5–100 mL


1–5 mL

Filter size


25 mm



13 mm

2. Pore size of filter membrane

For general purpose
HPLC columns packed with particles ≥ 3 µm, GC, SFC, ...


0.45 µm

Recommended for
UHPLC-, core-shell and HPLC columns, packed with particles ≤ 3 µm, GC, SFC, ...


0.20 µm

Product	Membrane pore size	Membrane dia.	Pack of	REF	List price per pack	Promotion price per pack *	Promotion price per filter *
Polytetrafluoroethylene (PTFE), for nonpolar liquids and gases							
CHROMAFIL® Xtra PTFE-20/13	0.20 µm	13 mm	100	729208	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra PTFE-20/25	0.20 µm	25 mm	100	729207	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra PTFE-45/13	0.45 µm	13 mm	100	729209	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra PTFE-45/25	0.45 µm	25 mm	100	729205	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra GF/PTFE-20/25	1.0/0.20 µm	25 mm	100	729270	~~$ 170.00~~	$ 60.00	$ 0.60
CHROMAFIL® Xtra GF/PTFE-45/25	1.0/0.45 µm	25 mm	100	729271	~~$ 170.00~~	$ 60.00	$ 0.60
Polyethersulfone (PES), for aqueous and slightly organic liquids							
CHROMAFIL® Xtra PES-20/25	0.20 µm	25 mm	100	729240	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra PES-45/25	0.45 µm	25 mm	100	729241	~~$ 135.00~~	$ 50.00	$ 0.50
Polyamide / Nylon (PA), for aqueous organic and aqueous medium polar liquids							
CHROMAFIL® Xtra PA-20/13	0.20 µm	13 mm	100	729248	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra PA-20/25	0.20 µm	25 mm	100	729212	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra PA-45/13	0.45 µm	13 mm	100	729249	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra PA-45/25	0.45 µm	25 mm	100	729213	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra GF/PA-20/25	1.0/0.20 µm	25 mm	100	729260	~~$ 170.00~~	$ 60.00	$ 0.60
CHROMAFIL® Xtra GF/PA-45/25	1.0/0.45 µm	25 mm	100	729261	~~$ 170.00~~	$ 60.00	$ 0.60
Cellulose acetate (CA), for filtration of water soluble oligomers and polymers, especially suited for biological macromolecules							
CHROMAFIL® Xtra CA-20/13	0.20 µm	13 mm	100	729254	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra CA-20/25	0.20 µm	25 mm	100	729226	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra CA-45/13	0.45 µm	13 mm	100	729255	~~$ 135.00~~	$ 50.00	$ 0.50
CHROMAFIL® Xtra CA-45/25	0.45 µm	25 mm	100	729227	~~$ 135.00~~	$ 50.00	$ 0.50

When ordering by e-mail or phone, please refer to

WIN26

Offer valid until December 31st, 2026

* Special promotion prices only valid for the mentioned CHROMAFIL® syringe filters starting at the first pack.



© kichigin19 – stock.adobe.com
Flyer Chromafil Xtra Winter Promotion en7a

www.mn-net.com

MACHEREY-NAGEL

MACHEREY-NAGEL Inc. · 924 Marcon Blvd., Suite 102 · Allentown, PA 18109 · USA

US +1 888 321 6224 sales-us@mn-net.com
DE +49 24 21 969-0 info@mn-net.com
CH +41 62 388 55 00 sales-ch@mn-net.com
FR +33 388 68 22 68 sales-fr@mn-net.com